1. American National Standard, Hemodialysis Systems (ANSI/AAMI RDS — 1992) (Association for the Advancement of Medical Instrumentation, Arlington, VA, 1993).

2. C. Sorber, W. Cooper, and E. Meyer, "Selection of a Field Method for Free Available Chlorine," *Disinfection-Water and Wastewater*, I. D. Johnson, Ed. (Ann Arbor Publishers, Ann Arbor, MI, 1975), pp. 91-112.

3. R. Bauer, B. E. Phillips and C. O. Rupp, "A Simple Test for Estimating Free Chlorine," *Journal AWWA* (November 1972), pp. 787-789.

4. J. Lieberman, N. M. Roscher, E. E. Meyer, and W. J. Cooper, "Development of the FACTS Procedure for Combined Forms of Chlorine and Ozone in Aqueous Solutions," *Environ Sci Technol* 14 (1980) 1395-1400.

5. "Amperometric Titration Method," *Standard Methods for the Examination of Water and Wastewater*, 18th Edition (American Public Health Association, Washington, D.C., 1991) 4-41 to 4-43.

**REFERENCES**

**HACH®  
STERICHEK®  
RESIDUAL CHLORINE  
REAGENT STRIPS**

**811900**

Hach Company  
100 Dayton Ave. • Ames, IA 50010 U.S.A.  
TOLL FREE: 888-ETS-STRIPS (1-888-387-7874)  
TEL: 970-278-4951 • FAX: 970-619-5025  
www.sterichек.com • etscustomerservice@hach.com

**DESCRIPTION**

SteriChek Residual Chlorine Reagent Strips provide a convenient, accurate means of measuring the concentration of chlorine bleach remaining in water being used to rinse out dialysate lines following disinfection of hemodialysis equipment.

A rapid screening method will detect levels above 0.5 ppm (mg/l) and a 30-second quantitative method allows estimation of concentrations between 0 and 5 ppm.

The qualitative rapid screening method can be used to determine that chlorine has been adequately rinsed from the machine. The quantitative method may be useful where corrective measures are undertaken on machines dispensing unacceptable levels of chlorine for extended periods and for testing containers, i.e. "bicarb jugs," disinfected with bleach.

Chlorine bleach detected in the rinse water is normally referred to as "residual chlorine." The results of the test indicate the level of chlorine bleach remaining in the dialysate lines at any point during the rinsing process.

The SteriChek Residual Chlorine Reagent Strips are supplied in a ready-to-use form. When placed in contact with the rinse water, according to directions, the test area changes color relative to the amount of chlorine remaining in the rinse. Results of 0.5 ppm or above indicate that further rinsing is necessary. <sup>1</sup>

**! WARNING**


- **Improper strip activation and color interpretation may result in patient injury.**
- **Keep all unused strips in the original bottle. Do not remove desiccant pack. Replace cap immediately and tightly after removing a strip; the strips must be protected from heat and humidity.**
- **Do not touch the reagent pad area. Do not allow the pad to come into contact with liquids or with work surfaces, as these may be contaminated with potentially interfering substances.**

**! IMPORTANT**


**Always compare test results to the color chart on the SteriChek bottle for proper interpretation.**

**Directions For Qualitative Results:**

1. Hold the reagent area in the rinse stream for 5 seconds.


2. If no color is apparent immediately upon removal, the chlorine level is less than 0.5 ppm and no further rinsing is indicated. Any pink/purple color indicates concentration of 0.5 ppm or greater.


**Directions For Quantitative Results:**

1. Fill standard sample cup with solution. Discard solution and re-fill.
2. Immerse reagent area of test strip in solution and move strip back and forth vigorously for 30 seconds.

3. Remove strip and compare to color chart within 10 seconds.

**To Retest a Sample:** Chlorine is consumed during the reaction. To retest a sample, empty the cup and repeat the entire procedure with a fresh sample.

**For Quality Control:**

Each facility should determine its own quality control procedure. Testing and recording test strip results with the control solution (see below) or with SteriChek Chlorine Control Tablets provides the user with a warning of a possible test strip error, potential use of outdated test strips, or of improperly stored or handled test strips. (SteriChek Chlorine Control Tablets are sold separately. Contact your distributor for more information.)

**Preparation of the Control Solution:**

Prepare the control solution by diluting concentrated chlorine bleach with Reverse Osmosis water.

**Dilution Guidelines:**

Stock Solution	Dilute 1 ml of chlorine bleach with 500 ml of Reverse Osmosis water. This stock solution contains approximately 100 ppm (mg/l) chlorine, has a useful life of 14 days and must be refrigerated.
Control Solution	Dilute 1 ml of the Stock Solution with 50 ml of Reverse Osmosis water. Use this control solution within 16 hours.

Measure the chlorine content of the control solution following the directions on the SteriChek product label. The assay value should fall within 1 to 3 ppm (mg/l).

Follow the procedure to prepare the stock solution with chlorine bleach typically used in your facility. Determine where your test strip or test kit values fall when using the control solution based on the typical chlorine bleach used in each facility. The assay value range of 1 to 3 ppm indicated here was obtained with concentrated chlorine solution that contained 6% sodium hypochlorite diluted to control solution level.

**Determination of Residual Iodine Levels**

The SteriChek Residual Chlorine Reagent Strips can be used to measure iodine in the absence of chlorine. For iodine concentration, follow the quantitative procedure and multiply the concentration level associated with the appropriate color block by 3.6.

**STORAGE**


The SteriChek Residual Chlorine Reagent Strips must be kept in the original bottle with the lid tightly closed to obtain the best results. Do not remove the desiccant pack. Store at temperatures between 60° - 90°F (16° - 32°C). Use within 6 months after first opening bottle. Do not use the test strips (from an opened or unopened bottle) after the expiration date.

**RESULTS**

The concentration of chlorine in rinse water is obtained by comparing the color of the reagent pad with color blocks on the label. The color blocks are calibrated in terms of chlorine concentration in parts per million (ppm = milligrams per liter or mg/l). Color blocks are designated as 0, 0.5, 1, 2 and 5 ppm (mg/l) chlorine. Concentrations (color development) which fall between color block values should be estimated.

**CHEMICAL PROPERTIES OF THE TEST**

SteriChek Residual Chlorine Reagent Strips react with both free chlorine and combined chlorine (monochloramines). Chlorine oxidizes syringaldazine (a colorless compound) to form a red-purple oxidation product. The reagent pad on the strip is buffered to pH 6.8 and contains potassium iodide. Combined chlorine in the form of monochloramines oxidizes the potassium iodide to iodine which in turn oxidizes syringaldazine to the colored form. <sup>2,3,4</sup>

**PERFORMANCE CHARACTERISTICS**

The performance characteristics of the SteriChek Residual Chlorine Reagent Strips are based on analytical studies using samples to which sodium hypochlorite was added to give a range of chlorine levels. Amperometric titration was used as the reference method for measuring chlorine levels. <sup>5</sup>

The sensitivity and accuracy of the strip test depends on several factors including variability of color perception, the variation in lighting condition and the possible presence of interfering substances. The lowest reliably detectable level of chlorine is about 0.25 ppm. Samples with reference chlorine concentrations falling between two color block values will give results ranging anywhere between those values. Results will generally be within one color block of the reference value. In blind studies, concentrations of 0.25 ppm have been read as less than 0.5 ppm, and concentrations of 1.0 ppm as greater than 0.5 ppm in all cases.

**LIMITATIONS**

The test will give a positive result with any substance which will oxidize syringaldazine directly or which will oxidize iodide to iodine under neutral pH conditions. These substances

include, among others, hypochlorite, chlorine, monochloramine, nitrogen trichloride, ozone, iodine, bromine and peroxide. SteriChek Residual Chlorine Reagent Strips are not suitable for testing for chlorine in bicarbonate concentrates.

Syringaldazine does not react with many of the substances found in non-purified water such as sulfate nitrate, chloride, copper, calcium and magnesium. Manganese (IV) and iron (III) give false positive results, while the presence of nitrate will yield a false negative. However, since rinsing of hemodialysis systems is usually done with highly purified water, it is unlikely that these potentially interfering substances will appear in the rinse water.

**AVAILABILITY**

Product Code 811900 SteriChek Residual Chlorine Reagent Test Strips includes five bottles of 100 reagent strips and a multilingual product manual. Also enclosed for your use are color-coded stickers that correspond to the color of the bottle label and kit box label. These stickers may be applied on the top of each bottle for easy product identification. Each sticker includes a space to record the date the bottle is opened.

These SteriChek testing products are also available from your distributor:

- 811902 0.1 ppm Total Chlorine DPD Kit
- 811903 0.1 ppm Total Chlorine DPD Refill Kit
- 811905 Residual Peroxide Reagent Strips
- 811906 Peracetic Acid Reagent Strips
- 811911 Sensitive 5 ppm Low-Range Hardness Strips
- 811912 Chlorine Control Tablets
- 811913 Residual Peroxide Control Tablets
- 811916 Bicarb pH Reagent Strips
- 812014 Blood Leak Reagent Strips

Made and Printed in the U.S.A. of US and imported content.

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