DESCRIPTION

The SteriChek Total Chlorine Test Kit provides a convenient, accurate means of measuring low levels of total chlorine (i.e. total chloramines plus free chlorine) in feed water used to prepare dialysate. The test indicates that the level of total chlorine is below the maximum allowable concentration of 0.1 ppm (mg/l) for chloramine. The chloramine level will always be less than or equal to the total chlorine level. ¹

The test kit will provide quantitative results in increments of 0.1 ppm utilizing the recognized standard method of analysis of the Association for the Advancement of Medical Instrumentation.² The kit includes reagent, apparatus and a quality control measure for monitoring results. A result of 0.1 ppm or greater indicates that the water tested should not be used to prepare dialysate and that carbon absorption media utilized to remove excess levels of chloramines may need to be replaced.²


⚠ WARNING

- Keep all unused powder pillows in the original package.
- Powder pillow contents may cause eye and respiratory tract irritation.
- Improper rinsing and storage of the sample tubes may impair the test results.
- **Do NOT throw tube caps away.**

DIRECTIONS

To detect low levels of total chlorine (chloramine).

1. Rinse the sample tubes thoroughly with water to be tested.
2. Fill each sample tube to the top line indicated on the tube with the water sample.
3. Open one DPD Total Chlorine Powder Pillow and add the contents to one of the sample tubes.
4. Cap the sample tube and shake gently for approximately 30 seconds.
5. Remove the caps from both tubes.
6. Place the sample tube with the prepared sample into the chamber on the inside of the comparator. Place the sample tube without the powder pillow in the outside chamber.
7. Hold the comparator to let light fall on the top of the tubes and rotate the color wheel inside the comparator to obtain a color match.
8. Read the total chlorine level as indicated on the comparator.
9. If the sample color falls between two colors, use the value of the higher color as the concentration.
10. IF THE COLOR WHEEL SHOULD BECOME WET, THE DISC COLORS MAY BECOME NON-UNIFORM. IF THAT SHOULD HAPPEN, REPLACE THE COMPARATOR.

For Quality Control:

Each facility should determine its own quality control procedure. Testing and recording results with the control solution (see below) or with SteriChek Chlorine Control Tablets provides the user with a warning of a possible operator error, potential use of outdated product, or of product improperly stored or handled. (SteriChek Chlorine Control Tablets are sold separately. Contact your distributor for more information.)

Preparation of the Control Solution:

Prepare the control solution by diluting concentrated chlorine bleach with Reverse Osmosis water.

Dilution Guidelines:

Stock Solution	Dilute 1 ml of chlorine bleach with 500 ml of Reverse Osmosis water. This stock solution contains approximately 100 ppm (mg/l) chlorine, has a useful life of 14 days and must be refrigerated.
Control Solution	Dilute 1 ml of the Stock Solution with 50 ml of Reverse Osmosis water. Use this control solution within 16 hours.

Measure the chlorine content of the control solution following the directions for the SteriChek Total Chlorine Test Kit. The assay value should provide a positive response within the range of the test kit.

Follow the procedure to prepare the stock solution with chlorine bleach typically used in your facility. Determine where your test kit values fall when using the control solution based on the typical chlorine bleach used in each facility. The positive response within the range of the test kit was obtained with concentrated chlorine solution that contained 6% sodium hypochlorite diluted to control solution level.

STORAGE

The SteriChek Total Chlorine powder pillows must be kept in their original packaging. Store at temperatures between 50° - 77°F (10° - 25° C). Do not use the powder pillows after the expiration date.

RESULTS

The concentration of total chlorine (total chloramines) in water is obtained by comparing the color of the prepared sample with the colors on the color wheel. The color wheel is calibrated in terms of total chlorine concentration in ppm (mg/l).

CHEMICAL PROPERTIES OF THE TEST

The SteriChek Total Chlorine Test Kit measures both free chlorine and combined chlorine (monochloramines and other chloramines). The reagent pillows contain DPD, potassium iodide (KI) and buffer. Chlorine oxidizes DPD (a colorless compound) to form a magenta (red) color. Chloramines oxidize potassium iodide to iodine which oxidizes DPD to form the magenta color. ¹

PERFORMANCE CHARACTERISTICS

The performance characteristics of the SteriChek Total Chlorine Test Kit are based on analytical studies using samples to which either sodium hypochlorite or monochloramine was added to give a range of free chlorine or combined chlorine levels.

LIMITATIONS

The test will give a positive result with any substance that will oxidize DPD directly or which will oxidize iodide to iodine at pH 6.6. These substances include, among others, chlorine, monochloramine, hypochlorite, nitrogen trichloride, ozone, iodine, bromine and peroxide.

The AAMI has published maximum allowable chemical contaminant levels for make-up water in hemodialysis systems.² These levels did not affect the results measured at 0.05, 0.1 and 0.2 mg/l monochloramine and 0.5, 0.6 and 0.7 mg/l free chlorine. Therefore, when rinsing of hemodialysis systems is done according to the AAMI Standard, no potentially interfering substances will appear in the rinse water. Oxidized manganese will give a false positive result.

AVAILABILITY

Product Code 811902 SteriChek Total Chlorine Test Kit includes one bag of 100 DPD powder pillows, two sample tubes, one Comparator including color wheel, and a multilingual product manual.
Product Code 811903 SteriChek Total Chlorine Test Kit Refill includes five bags of 100 DPD powder pillows and a multilingual product manual.

These SteriChek testing products are also available from your distributor:

- 811900 Residual Chlorine Reagent Strips
- 811905 Residual Peroxide Reagent Strips
- 811906 Peracetic Acid Reagent Strips
- 811911 Sensitive 5 ppm Low-Range Hardness Strips
- 811912 Chlorine Control Tablets
- 811913 Residual Peroxide Control Tablets
- 811916 Bicarb pH Reagent Strips
- 812014 Blood Leak Reagent Strips

Made and Printed in the U.S.A. of US and imported content.

REFERENCES

1. C. Sorber, W. Cooper and G. Water, "Selection of a Field Method for Free Available Chlorine," Disinfection-Water and Wastewater, L.D. Johnson, Ed. (Ann Arbor publishers, Ann Arbor, MI, 1975), p.91-112.
2. 1993 Association for the Advancement of Medical Instrumentation (ANSI/AAMI:RD5-1992) p.5.
3. "Amperometric Titration Method", Standard methods for the Examination of Water and Wastewater, 18th Edition (American Public Health Association, Washington, DC, 1991) p.4-41 to 4-43.

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